

## BATEC COMPATIBILITY CHART 2015

---

### 1. INTRODUCTION

---

With our BATEC DOCKING-BAR (Pat. Pending) all our range of Batecs can be installed in most of the wheelchairs –rigid or folding frame, folding frame with removable footplate. Thanks to our set of clamps and brackets a Batec can be fitted to different sorts of frames and a tubes sizes.

In custom-made wheelchairs the installation might be different in every specific case.

You will now find some information about:

- Minimum and maximum sizes that a wheelchair should have to install a Batec.
- A chart of the most common wheelchairs in which a Batec can be fitted including the appropriate clamps and brackets.
- A chart of non-compatible wheelchairs.
- A chart of wheelchair manufacturers who are welding our BATEC DOCKING-BAR to some of their wheelchairs.

**! IF YOUR WHEELCHAIR DOES NOT APPEAR IN THIS LIST, PLEASE CONTACT US AND WE WILL LOOK FOR THE BEST SOLUTION.**

### 2. MINIMUM AND MAXIMUM MEASURES

---

SEAT WIDTH	30 to 50 cm
MINIMUM HEIGHT FOOTREST – FLOOR	5 cm
MINIMUM FRONT SEAT HEIGHT	45 cm
MINIMUM DISTANCE SEAT - FOOTREST	30 cm
MAXIMUM USER WEIGHT	110 kg

### 3. COMPATIBLE WHEELCHAIRS

BRAND	MODEL	BRACKETS	CLAMPS
COLOURS	SPAZZ / SPAZZ G	RIGID	D2 25
COLOURS	SABER	RIGID	D2 25
COLOURS	SHOCKBLADE	RIGID	D2 25
COLOURS	RAZORBLADE	RIGID	D2 25
COLOURS	BOING!	RIGID	D2 25
COLOURS	ZEPHYR	RIGID	D2 19
COLOURS	ECLIPSE	RIGID	D2 25
COLOURS	GENESIS	FOLDABLE	D2 30
CYCLONE	XLR-8	RIGID	D2 C
CYCLONE	OTHER MODELS	RIGID	! The installation depends on the configuration of every particular unit
GTM	GTM1	RIGID	D2 25
GTM	BASIC	RIGID	D2 25
GTM	HAMMER	RIGID	D2 35
GTM	CARISA	RIGID	D2 19
GTM	CHALLENGER	RIGID	D2 35
GTM	SHOCK	RIGID	D2 35
GTM	PRESIDENT	RIGID	D2 19
GTM	ENDEAVOUR	RIGID	D2 35
GTM	MUSTANG	RIGID	D2 30
KUSCHALL	K4	RIGID	D2 30
KUSCHALL	K4 AIRLITE	RIGID	D2 30
KUSCHALL	KSL	RIGID	D2 30
KUSCHALL	R33	RIGID 50°	D2 30
KUSCHALL	ADVANCE	RIGID	D2 A
KUSCHALL	COMPACT	FOLDABLE	D2 25
KUSCHALL	ULTRALIGHT	FOLDABLE	D2 25
KUSCHALL	CHAMPION	FOLDABLE	D2 30
MEYRA	HURRICANE	RIGID	D2 25
MEYRA	XSTAR	RIGID	D2 25
MEYRA	X1	FOLDABLE	D2 25
MEYRA	ZX1	FOLDABLE	D2 25
MEYRA	ZX3	FOLDABLE	D2 25
MEYRA	AVANTI	FOLDABLE	D2 25
ORACING	DIFFERENT MODELS	RIGID	! The installation depends on the configuration of every particular unit
OTTOBOCK	VENTUS	RIGID	D2 30
OTTOBOCK	BLIZZARD	RIGID	D2 B
OTTOBOCK	MOTUS CV	FOLDABLE	D2 25
OTTOBOCK	MOTUS CS	FOLDABLE	D2 25
OTTOBOCK	AVANTGARDE	FOLDABLE	D2 25

PANTHERA	S2	RIGID	D2 19
PANTHERA	U2	RIGID	D2 19
PANTHERA	U2 LIGHT	RIGID 50°	D2 19
PROGEO	JOKER	RIGID	D2 J
PROGEO	JOKER ENERGY	RIGID	D2 J
PROGEO	JOKER EVOLUTION	RIGID	D2 J
PROGEO	NOIR	RIGID	D2 T
PROGEO	TEKNA	FOLDABLE	D2 T
PROGEO	TEKNA ADVANCE	FOLDABLE	D2 T
PROGEO	BASIC LIGHT	FOLDABLE	D2 22
PROGEO	EXELLE VARIO	FOLDABLE	D2 22
PROGEO	EXELLE	FOLDABLE	D2 22
QUICKIE	HELIUM	RIGID	D2 H
QUICKIE	HELIUM PRO	RIGID	D2 H
QUICKIE	TI	RIGID	D2 25
QUICKIE	ARGON	RIGID	D2 30
QUICKIE	ARGON TI	RIGID	D2 30
QUICKIE	REVOLUTION	RIGID	D2 25
QUICKIE	LIFE R	RIGID	D2 22
QUICKIE	XENON	FOLDABLE	D2 H
QUICKIE	NEON	FOLDABLE	D2 30
QUICKIE	NEON REP. EXTR.	FOLDABLE	D2 25
QUICKIE	LIFE	FOLDABLE	D2 22
QUICKIE	EASY MAX	FOLDABLE	D2 22
QUICKIE	HELIX	FOLDABLE	D2 22
RGK	DIFFERENT MODELS	RIGID	! The installation depends on the configuration of every particular unit
TILITE	Z10	RIGID	D2 30
TILITE	ZRA	RIGID	D2 30
TILITE	ZR	RIGID	D2 30
TILITE	TR	RIGID	D2 30
TILITE	AEROZ	RIGID	D2 30
TILITE	AEROR	RIGID	D2 25
TILITE	AEROX	FOLDABLE	D2 25
TILITE	TX	RIGID	D2 25
TILITE	2GX	RIGID	D2 25
TOP-END	TERMINATOR STD/TITAN.	RIGID	D2 25
TOP-END	CROSSFIRE	RIGID	D2 25
WOLTURNUS	DIFFERENT MODELS	RIGID AND FOLDABLE	! The installation depends on the configuration of every particular unit

#### 4. NON COMPATIBLE WHEELCHAIRS

---

BRAND	MODEL	REASON OF INCOMPATIBILITY
KÜSCHALL	K4 CARBONO	Frame material
KÜSCHALL	CHAMPION CARBONO	Frame material
MEYRA	XR	Design of the frame
PANTHERA	X	Frame material
PROGEO	YOGA	Design of the folding mechanism
PROGEO	EGO	Design of the folding mechanism
QUICKIE	XENON REP EXTR	Design of the folding mechanism
TILITE	AERO T	Design of the frame

#### 5. WHEELCHAIR MANUFACTURERS AND BATEC

---

Some of best wheelchair manufacturers agreed to weld the BATEC DOCKING-BAR to their wheelchairs. Thanks to that, a perfect commitment between strength, lightness and aesthetics is achieved.

Nowadays, our welded BATEC DOCKING-BAR can be found in the following wheelchairs:

**RGK:**

Units 8a/b/c. Ring Road, Zone 2  
Burntwood Business Park. Staffordshire  
WS7 3JQ  
United Kingdom  
T: (+49) 1543670077  
<http://www.rgklife.com/en/>  
[info@RGKlife.com](mailto:info@RGKlife.com)

**GTM-MOBIL MAREK-TECHNOLOGY:**

WALISIEWICZ SP.J  
UL.DORODNA 16  
03-195 WARSZAWA  
T: (+48) 508364720  
[www.gtmmobil.com.pl](http://www.gtmmobil.com.pl)  
[info@gtmmobil.com.pl](mailto:info@gtmmobil.com.pl)

**WOLTURNUS:**

Halkærvej 24B  
DK-9240 Nibe  
Denmark  
T: (+45) 96717170  
<http://www.wolturnus.com/>  
[info@wolturnus.dk](mailto:info@wolturnus.dk)

**CYCLONE MOBILITY:**

Orchard House,  
Sunk Island Road,  
Ottringham,  
Hull, East Yorkshire,  
HU12 0DX  
T: (+49) 8001804850  
[www.cyclonemobility.com](http://www.cyclonemobility.com)  
[info@cyclonemobility.com](mailto:info@cyclonemobility.com)

**! OTHER WHEELCHAIR MANUFACTURERS ARE BEING ABLE TO WELD THE BATEC DOCKING BAR ON OTHER MODELS IN A FUTURE. ASK FOR IT WHEN ORDERING YOUR NEW WHEELCHAIR.**